

Written By: Robert Stillwell

Written For: NCAR

This file is used to define the initial state of the water vapor DIAL Main labview program when idling with no children loaded

These are the actual names that will appear in the tabs for the container.

The tag "Unused" tells the container not to populate a tab. Note that this

command is case sensitive.

Tab Names;;

Main Controls; Unused; Unused; Unused; Unused; Unused; Unused; Unused; Unused; Unused; Unused; Unused;;

In its current state, these two variables are the width and height of the

tab containers. This location is, however, a catch all for data you need to

read in quickly from this configure file

Global Variables;;

540;985;;

This is the relative file path(s) (identified by the .\) to the program(s)

that the container should populate the tabs with

Relative File Paths;;

;;

These numbers correspond to the TypeDef_MainControlSystem.ctl. They tell

the main container what to populate the tabs with.

Function Types;;

;;

the paths to logging for main container

Logging Paths;;

.\Data\Container\;

ContainerLogging;;

the name of the dial unit used to grab the correct config files

Unit Name;;

DIAL2;;